


# ArterioGlide

Ref V-AC20G4.5



Arterial Catheterization Kit

# Arterial Blood Pressure and blood gas analysis


- ▶ **Guidewire:** made of nitinol, with characteristics, unique superelastic, antikinking and shape memory properties.

- ▶ **Great biocompatible material:** catheter made of PU (polyurethane)

- best surface lubricity
- Flexibility antikinking
- Antistick properties


- ▶ **Radiopaque:** this facilitates quick visualization under X-ray.

- ▶ **Soft and tapered tip:** this design facilitates smooth insertion, thin-delicate needle and reduces vessel trauma during insertion.




# ArterioGlide

**Assures** the medical staff an easy, safe and atraumatic punctures of the peripheric arteries, the Insertion is performed using **the seldinger technique** and the kit has the necessary accesories for the performance of the procedure.



**The intention of** use of the **ArterioGlide** Arterial Catheter is measurement of blood pressure and blood gas analysis, and blood collection avoiding nonconformities in the patient.

**ArterioGlide** can be used in the following arteries such as: radial artery, brachial artery, femoral artery and pedis artery.



Part Nro	Introducer Needle	Arterial Catheter	Guide Wire Diameter
V-AC20G	20 Ga x 1-1/2" (4cm)	20 Ga x 1-3/4" (4.5cm)	0.21" x 1-3/4" (4,5cm)



**HOSMED-INC**  
**3403 NW 82nd Avenue. Suite #102.**  
**Doral, Florida 33122. USA**  
**[www.hosmed-inc.com](http://www.hosmed-inc.com)**  
**P: (954) 505-1432**  
**E-mail: [marketing@hosmed-inc.com](mailto:marketing@hosmed-inc.com)**

**H O S M E D**