



Hydraulic Stretcher NS5H



Features



Hydraulic treadle for rise



(1 Leg section lowers down to make head higher (Than Leg Section)
 2. Back section lowers down to make leg section higher (Than head Section)
 3. Step on the middle of the pedal, make the bed base fall down steadily



Four double side casters of 6 inch diameter, simultaneously lockable by the step pedal at both sides.

TPR tire no worn out after running 30KM, save anti-winding hard shell, united forming without bolts



Waterproof PU mattress, convenient for transfer patient



Foldable aluminium alloy guardrail



Bumper wheel detachable, max. protect the bed and patient during transportation



Full SS IV pole, bearing 15kg

Technical Chart

QUALITY CONTROL





Technical parameters

| | |
|-----------------------|----------------------|
| Size | 2100*850*(630-930)mm |
| Back-rest | 0-75° (±10°) |
| Knee-rest | 0-45° |
| Reverse Trendelenburg | 0-15° |
| Trendelenburg | 0-15° |

Technical configuration

| | |
|------------------------------------|--------|
| USA Hydraulic Pump | 2 pcs. |
| Crank System | 1 set |
| Back Rest Air Spring | 2 pcs. |
| 8" German Castors, Central Locking | 4 pcs. |
| Emergency Push Handle | 1 set |
| Folding Side Rail | 1 set |
| 74mm Wall Bumper Wheel | 1 set |
| IV pole prevision | 2 pcs. |
| Mattress | 1 pcs. |
| IV Pole | 1 set |
| Central Locking Break Pedal | 1 pcs. |
| X-Ray Translucent Bed Platform | 1 set |





3403 NW 82nd Avenue, Suite 102 E.
Doral, Florida 33122. USA
www.hosmed-inc.com