


Rover Transradial Introducer Kit



* Longer taper ensure smooth transition between sheath and dilator

* Excellent pushability, flexibility and radiopacity

Cod	Type	Sheath ID(mm)L(mm)		Dilator ID(mm)L(mm)		Guide wire OD	Needle (G/mm)	Hydrophilic Guide wire	Hydrophilic Coating
V-RRO5F2I	5F	1.86	110	1.1	170	0.021"	21/45	No	No
V-RRO6F2I	6F	2.18	110	1.1	170	0.021"	21/45	No	No



3403 SW 82ND AVENUE. SUITE #102.

DORAL, FLORIDA 33122. USA

WWW.HOSMED-INC.COM

P: +1 786 464 0361

E-mail: marketing@hosmed-inc.com